

# MicroVal and NordVal certification of Ready-To-Use Culture Media, Easy Plate EC for Enumeration of *E. coli* and Coliforms

## in a Broad Range of Foods, Pet Foods and Environmental Samples

Daiki Utsunomiya<sup>1</sup>, Suzanne J. Jordan<sup>2</sup>, Kentaro Takenaka<sup>3</sup>, Mai Shimizu<sup>3</sup>, Takeo Suzuki<sup>3</sup> and Shinichiro Sugiura<sup>1</sup>

(1)Kikkoman Biochemifa Company, (2)Campden BRI, (3)Kikkoman Corporation




### Introduction

- **Easy Plate** series (AC, EC, EC, EB, SA and YM-R) are Ready-to-Use (RTU) medium which are manufactured and sold by Kikkoman Biochemifa Company. Compared to conventional plate media, Easy Plate provides various advantages such as **time reduction, ease of operation, space-saving and reduction of plastic.**

- **Easy Plate EC (E-EC)** can be used for *E. coli* and coliforms.

*E. coli* : blue colony  
Coliforms : purple colony



- According to ISO (16649-2:2001 and 4832:2006), tryptone bile X-glucuronic agar medium (TBX) and violet red bile lactose agar medium (VRBL) are used as conventional methods for *E. coli* and coliforms, respectively.

- In this study, the experiments were conducted to evaluate the specificity, selectivity, repeatability, accuracy and relative trueness of the E-EC for the enumeration of *E. coli* and coliforms as required by ISO 16140-2:2016.

### Methods

The method validation study was done according to ISO 16140-2:2016, using ISO 16649-2:2001 and ISO 4832:2006 as the reference methods. Analysis with E-EC was performed following manufacturer's instructions.

#### 1. Inclusivity and Exclusivity


Pure cultures of known provenance were tested for the inclusivity and the exclusivity. Each test was performed once with E-EC, TBX (for *E. coli*), VRBL (for coliforms) and non-selective agar.

#### 2. Accuracy profile

Seven categories (Milk and dairy products, Fishery products, Produce and fruits, Multi-components foods, Raw and RTC meat and poultry, Pet food and animal feed and Environmental samples) were tested for accuracy profile study.

#### 3. Relative trueness

Relative trueness study was conducted using a combination naturally and artificially contaminated samples across seven categories (same as the accuracy profile).



### Result 1. Inclusivity and Exclusivity

For the inclusivity study, a total of 50 pure cultures of *E. coli* and 51 of coliforms were tested. Forty-nine of the 50 *E. coli* and 50 of 51 coliform strains gave typical *E. coli* colonies on E-EC. For the exclusivity study, a total of 30 pure non-target cultures were tested separately for E-EC versus TBX and E-EC versus VRBL. E-EC and TBX showed similar selectivity for *E. coli*. Otherwise, E-EC showed better selectivity compared to VRBL for coliforms (False-positive strains are 1 and 6, respectively). **These results indicated that E-EC has similar specificity for *E. coli* as TBX and had a greater selectivity for coliforms than VRBL.**

**Table 1.** Summary of inclusivity and exclusivity

	E. coli	Coliforms	Non target	False positive
E-EC	49/50	50/51	57/60	<i>Shigella sp.</i>
ISO methods	50/50	51/51	52/60	<i>Shigella sp.</i> <i>Vibrio sp.</i>

**Table 2.** Result of the inclusivity study

Strains	No. of Strains	E-CC*	TBX*	VRBLA*
<i>E. coli</i>	50	49	50	-
<i>Citrobacter</i>	10	9	-	10
<i>Enterobacter</i>	19	19	-	19
<i>Escherichia</i>	5	5	-	5
<i>Klebsiella</i>	6	6	-	6
<i>Cronobacter</i>	4	4	-	4
<i>Hafnia</i>	1	1	-	1
<i>Kluyvera</i>	1	1	-	1
<i>Serratia</i>	2	2	-	2
<i>Siccibacter</i>	1	1	-	1
<i>Franconibacter</i>	2	2	-	2

\* Number of strains giving anticipated results

**Table 3.** Result of the exclusivity study

No.	Organisms	Source or Identity	E-EC	VRBL	BGLB
1	<i>Aeromonas salmonicida</i>	NCTC 10402	✓	✓	NT
2	<i>Bacillus cereus</i>	ATCC 10876	✓	✓	NT
3	<i>Bacillus subtilis</i>	ATCC 6633	✓	✓	NT
4	<i>Carnobacterium divergens</i>	Brie	✓	✓	NT
5	<i>Edwardsiella tarda</i>	NCTC 10396	✓	✓	NT
6	<i>Flavobacterium species</i>	Bamboo shoots	✓	✓	NT
7	<i>Lactobacillus acidophilus</i>	industrial	✓	✓	NT
8	<i>Listeria monocytogenes</i>	Soft cheese	✓	✓	NT
9	<i>Aeromonas hydrophila</i>	tin of milk	✓	✓	NT
10	<i>Pasteurella bettyae</i>	NCTC 10535	✓	✓	NT
11	<i>Pediococcus pentasaceus</i>	Brine	✓	✓	NT
12	<i>Proteus mirabilis</i>	Poultry	✓	✓	NT
13	<i>Proteus vulgaris</i>	Poultry	✓	✓	NT
14	<i>Providencia rettgeri</i>	NCTC 7475	✓	✓	NT
15	<i>Pseudomonas aeruginosa</i>	NCIMB 10753	✓	✓	NT
16	<i>Pseudomonas fluorescens</i>	NCIMB 10586	✓	✓	NT
17	<i>Pseudomonas fluorescens</i>	water	✓	✓	NT
18	<i>Salmonella Typhimurium</i>	ATCC 14028	✓	✓	NT
19	<i>Acinetobacter calcoaceticus</i>	Sesame seeds	✓	✓	NT
20	<i>Acinetobacter lwoffii</i>	Tomatoes	✓	✓	NT
21	<i>Burkholderia gladioli</i>	Industrial	✓	✓	NT
22	<i>Shewanella putrefaciens</i>	NCTC 10736	✓	✓	NT
23	<i>Shigella boydii</i>	NCTC 11312	✓	✗	✓
24	<i>Shigella flexneri</i>	NCTC 9950	✓	✗	✓
25	<i>Shigella sonnei</i>	ATCC 25931	✗	✗	✓
26	<i>Staphylococcus aureus</i>	NCIMB 12702	✓	✓	NT
27	<i>Vibrio mimicus</i>	NCTC 11435	✓	✗	✓
28	<i>Vibrio parahaemolyticus</i>	NCTC 11344	✓	✗	✓
29	<i>Yersinia enterocolitica</i>	NCTC 10460	✓	✗	✓
30	<i>Morganella morganii</i>	Mince	✓	✓	NT

Total 29/30 24/30


“✓” indicates successfully suppressed  
“✗” indicates false-positive. NT: Not tested

### Result 2. Accuracy profile


A total of 210 samples (comprising two items per category at three different contamination levels) were used in this study. The accuracy of the E-EC was satisfied as all 7 categories met the 0.5log AL or the re-calculated AL for both *E. coli* and coliforms.

### Result 3. Relative trueness

A total of 105 samples across seven categories were tested for the relative trueness study. **The results showed that no significant difference between E-EC and the ISO methods for all seven categories.** For *E. coli*, the correlation coefficient between E-EC and TBX for all categories was 0.984 as shown in Figure 3-1. For coliforms, the correlation coefficient between E-EC and VRBL was 0.979 as shown in Figure 3-2.



**Figure 3-1.** Scatter plot of the reference method versus E-EC results for *E. coli*.



**Figure 3-2.** Scatter plot of the reference method versus E-EC results for coliforms.

### Conclusion

- Easy Plate EC showed similar inclusivity to the ISO methods (TBX and VRBL), and in term of exclusivity, Easy Plate EC differentiated more non-target organisms compared to VRBL.
- Easy Plate EC showed satisfactory results for accuracy profile and high correlation with TBX and VRBL.
- These results suggest that **Easy Plate EC is a comparative method to the ISO reference methods and applicable to a broad range of samples** for enumeration of both *E. coli* and coliforms.

### Contact information

Daiki Utsunomiya

✉ [dutsunomiya@mail.kikkoman.co.jp](mailto:dutsunomiya@mail.kikkoman.co.jp)

Kikkoman Biochemifa Company

<https://biochemifa.kikkoman.com/e/kit/>

LinkedIn



LinkedIn



Website

